

Kód	ks v balení	obor	názov	vlákno						ihla				špecifikácia	trieda
				Monofilament / Braid	materiál vlákna	farba	USP	dec	dĺžka (cm)	SA/DA	zakrivenie	ihla	mm		
12D10A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	5/0	1	45	Single Armed	3/8	ExtraCut (XT)	11mm	Low-trauma reverse cutting needle	III
12N15G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	1/2	Taper	22mm	ExtraBlack Needle	III
12N20G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	1/2	Taper	22mm	ExtraBlack Needle	III
12S07A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	6/0	0.7	75	Single Armed	3/8	Reverse Cutting	11 mm		III
12S10B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	5/0	1	75	Single Armed	3/8	Reverse Cutting	13 mm		III
12S10C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	5/0	1	75	Single Armed	3/8	Tapercutting (KL)	14 mm		III
12S10D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	5/0	1	75	Single Armed	3/8	Reverse Cutting	16 mm		III
12S10G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	5/0	1	75	Single Armed	3/8	Reverse Cutting	19 mm		III
12S15A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	Straight	Reverse Cutting	50 mm		III
12S15B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	24 mm		III
12S15C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Tapercutting (KL)	18 mm		III
12S15D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	13 mm		III
12S15E	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	16 mm		III
12S15F	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	19 mm		III
12S15G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	1/2	Taper	22mm		III
12S15H	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Tapercutting (KL)	16 mm		III
12S15R	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	3/8	Tapercutting (KL)	16 mm		III
12S15T	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	4/0	1.5	75	Single Armed	1/2	Taper	16 mm		III
12S20A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	Straight	Reverse Cutting	50 mm		III
12S20B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	3/8	Reverse Cutting	24 mm		III
12S20C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	3/8	Tapercutting (KL)	18 mm		III
12S20D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	90	Single Armed	3/8	Reverse Cutting	30 mm		III
12S20E	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	3/8	Reverse Cutting	16 mm		III
12S20F	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	3/8	Reverse Cutting	19 mm		III
12S20G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	1/2	Taper	22mm		III
12S20K	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	1/2	Reverse Cutting	22mm		III
12S20L	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	1/2	Taper	30 mm		III
12S20Q	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	3/0	2	75	Single Armed	1/2	Tapercutting (KL)	22mm		III
12S30A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	90	Single Armed	1/2	Tapercutting (KL)	36 mm	Fine Needle	III
12S30AL	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	90	Single Armed	1/2	Tapercutting (KL)	36 mm		III
12S30B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	Straight	Reverse Cutting	50 mm		III
12S30C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	3/8	Reverse Cutting	24 mm		III
12S30D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	90	Single Armed	3/8	Reverse Cutting	30 mm		III
12S30E	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	3/8	Reverse Cutting	19 mm		III
12S30F	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	90	Single Armed	1/2	Reverse Cutting	36 mm		III
12S30G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	1/2	Reverse Cutting	26 mm		III
12S30J	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	1/2	Taper	26 mm		III
12S30N	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	2/0	3	75	Single Armed	1/2	Taper	30 mm		III
12S35A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	1/2	Tapercutting (KL)	36 mm		III
12S35B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	75	Single Armed	Straight	Reverse Cutting	50 mm		III
12S35C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	3/8	Reverse Cutting	30 mm		III
12S35D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	75	Single Armed	1/2	Taper	36 mm		III

12S35E	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	1/2	Tapercutting (KL)	48 mm	III
12S35F	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	1/2	Taper	40 mm	III
12S35G	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	1/2	Reverse Cutting	36 mm	III
12S35H	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	0	3.5	90	Single Armed	1/2	Taper	48 mm	III
12S40A	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Reverse Cutting	40 mm	III
12S40B	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Tapercutting (KL)	40 mm	III
12S40C	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Taper	40 mm	III
12S40D	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Taper	48 mm	III
12S40F	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	3/8	Reverse Cutting	36 mm	III
12S40H	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Reverse Cutting	48 mm	III
12S40N	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Tapercutting (KL)	36 mm	III
12S40P	36	Absorbable	OPTIME R	Braid	PGA Rapid	Undyed	1	4	90	Single Armed	1/2	Tapercutting (KL)	45 mm	III