

Kód	ks v balení	obor	názov	vlákno						ihla				špecifikácia	trieda
				Monofilament / Braid	materiál vlákna	farba	USP	dec	dĺžka (cm)	SA/DA	zakrivenie	ihla	mm		
91120A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	White	3/0	2	75	Single Armed	3/8	Taper	18 mm	III	
91130C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	White	2/0	3	75	Single Armed	1/2	Taper	26 mm	III	
91170C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	White	5	7	75	Single Armed	1/2	Tapercutting (KL)	50 mm	III	
91R15A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	4/0	1.5	250	reel		0		III	
91R20A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	250	reel		0		III	
91R30A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	250	reel		0		III	
91R35A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	250	reel		0		III	
91R40A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	250	reel		0		III	
91R50A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	250	reel		0		III	
91R60A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	250	reel		0		III	
91R70A	24	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	5	7	250	reel		0		III	
91S10A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	5/0	1	75	Single Armed	3/8	Reverse Cutting	16 mm	III	
91S15A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	12 mm	III	
91S15D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	4/0	1.5	75	Single Armed	3/8	Reverse Cutting	20 mm	III	
91S15E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	4/0	1.5	75	Single Armed	Straight	Reverse Cutting	50 mm	III	
91S20B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	Straight	Reverse Cutting	50 mm	III	
91S20C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	1/2	Reverse Cutting	30 mm	III	
91S20D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	1/2	Taper	30 mm	III	
91S20E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	1/2	Reverse Cutting	22mm	III	
91S20H	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	3/8	Reverse Cutting	19 mm	III	
91S20J	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	3/8	Taper	16 mm	III	
91S20K	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	3/8	Reverse Cutting	24 mm	III	
91S20L	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	1/2	Taper	26 mm	III	
91S20N	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3/0	2	75	Single Armed	Straight	Reverse Cutting	60 mm	III	
91S30A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	3/8	Reverse Cutting	19 mm	III	
91S30B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	3/8	Reverse Cutting	26 mm	III	
91S30C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	1/2	Taper	26 mm	III	
91S30D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	1/2	Reverse Cutting	30 mm	III	
91S30E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	Straight	Reverse Cutting	50 mm	III	
91S30F	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	1/2	Reverse Cutting	36 mm	III	
91S30G	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	1/2	Taper	30 mm	III	
91S30H	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	1/2	Reverse Cutting	26 mm	III	
91S30K	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	3/8	Taper	18 mm	III	
91S30L	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	Straight	Reverse Cutting	36 mm	III	
91S30M	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	3/8	Reverse Cutting	30 mm	III	
91S30N	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2/0	3	75	Single Armed	Straight	Reverse Cutting	60 mm	III	
91S35A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Reverse Cutting (OHS)	22mm	OrthoShape extra strong needle III	
91S35B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	3/8	Reverse Cutting	26 mm	III	
91S35C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Reverse Cutting	30 mm	III	
91S35D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	3/8	Reverse Cutting	30 mm	III	
91S35E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Reverse Cutting	36 mm	III	
91S35F	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Taper	26 mm	III	
91S35G	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	Straight	Reverse Cutting	50 mm	III	
91S35J	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Reverse Cutting	26 mm	III	
91S35K	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Taper	30 mm	III	

91S35M	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Taper	36 mm		III
91S35P	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Tapercutting (KL)	26 mm		III
91S35Q	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	3/8	Taper	30 mm		III
91S35R	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	0	3.5	75	Single Armed	1/2	Tapercutting (KL)	48 mm		III
91S40A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting (OHS)	22mm	OrthoShape extra strong needle	III
91S40B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting	26 mm		III
91S40C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting	30 mm		III
91S40D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Reverse Cutting	30 mm		III
91S40E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting	36 mm		III
91S40F	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting	40 mm		III
91S40G	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Reverse Cutting	48 mm		III
91S40H	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Taper	26 mm		III
91S40J	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Reverse Cutting	26 mm		III
91S40K	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	Straight	Reverse Cutting	50 mm		III
91S40L	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Taper	36 mm		III
91S40N	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	1/2	Taper	30 mm		III
91S40R	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Taper	40 mm		III
91S40T	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Reverse Cutting	36 mm		III
91S40U	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Reverse Cutting	40 mm		III
91S40V	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	Straight	Reverse Cutting	60 mm		III
91S40W	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	1	4	75	Single Armed	3/8	Taper	30 mm		III
91S50A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Reverse Cutting (OHS)	22mm	OrthoShape extra strong needle	III
91S50B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Reverse Cutting	30 mm		III
91S50C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	3/8	Reverse Cutting	30 mm		III
91S50E	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Reverse Cutting (OHS)	40 mm	OrthoShape extra strong needle	III
91S50J	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Tapercutting (KL)	48 mm		III
91S50K	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Taper	40 mm		III
91S50L	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Reverse Cutting	36 mm		III
91S50M	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	2	5	75	Single Armed	1/2	Reverse Cutting	40 mm		III
91S60A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Reverse Cutting	26 mm		III
91S60B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Reverse Cutting (OHS)	36 mm	OrthoShape extra strong needle	III
91S60C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Reverse Cutting	40 mm		III
91S60D	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Reverse Cutting	48 mm		III
91S60G	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	3/8	Reverse Cutting	36 mm		III
91S60H	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Reverse Cutting (OHS)	36 mm	OrthoShape extra strong needle	III
91S60J	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	3	6	75	Single Armed	1/2	Taper	36 mm		III
91S70A	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	5	7	75	Single Armed	1/2	Reverse Cutting (OHS)	40 mm	OrthoShape extra strong needle	III
91S70B	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	5	7	75	Single Armed	3/8	Reverse Cutting	50 mm		III
91S70C	36	Others non absorbable	POLYTRESSE	Braid	POLYESTER	Green	5	7	75	Single Armed	1/2	Tapercutting (KL)	50 mm		III